- 18)Don't use the device in the environment with higher than 40% or lower than 5%, and over RH 15%- 90%.
19)The protective glass outside the LCD frame is very important and also the fragile part of the instrument, so it must be used carefully.
20)Please do not charge the non-rechargeable battery and do not throw the battery into the fire.
21)Please do not expose the product to the sun or touch water.
22)To help avoid strangulation, keep the air tube away from infants, toddlers or children.
23)Keep components out of the reach of infants, toddlers or children. This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers or children.
24)Keep batteries out of the reach of infants, toddlers or children. If the monitor and other components are wet:
- In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
- In locations exposed to vibrations or shocks.

3.Potential adverse reaction

- 1)Patients with severe blood circulation disorders and blood diseases, inflate the arm cuff until necessary may cause blood stasis or paralysis of the arm;
2)Patients with severe blood circulation disorders and blood diseases may cause acute internal bleeding due to the compression of the arm.

4.Waste disposal

- 1)Direct disposal of electronic products and batteries in garbage cans will cause harm to the environment. Please dispose them according to the local laws.
2)Don't throw the Electronic Sphygmomanometers in the garbage can at the end of use; please dispose it according to the local laws or contact with the manufacturer for recycling.

6

- 3)Wrap the arm Cuff firmly in place around your left upper arm. It is advisable to tie the Cuff properly that the arm will not feel physical tension.

- The bottom edge of the arm cuff should be 2 to 3 cm above the elbow. Air tube is on the inside of your arm and aligned with your middle finger.

- 4)Secure closed with the fabric fastener.
Notes:
• When you take a measurement on the right arm, the

Table of Contents

I. Preface.....	.1
II. Range of application.....	.2
III. Important Safety Information.....	.3
IV. Features and Specification.....	.7
V. Structure and descriptions of product.....	.8
VI. Settings.....	.10
VII. Preparation.....	.12
VIII. Method to measure the blood pressure.....	.16
IX. Requirements of Storage and maintenance.....	.19
X. Calibration and Service.....	.20
XI. Instructions for Waste Disposal.....	.20
XII. Symbol Description.....	.22
XIII. Service information.....	.23
XIV. EMC Declaration.....	.23
XV. Manufacture.....	.28
XVI. European Representative.....	.28

IV.Features and Specification


1.List of product

No.	Part Name	Model	Quantity
1)	Main unit	X5	1 set
2)	Cuff	URN-224Z	1 set
3)	Instructions Manual	/	1 copy

2.Technical specification

Device Name	Electronic Sphygmomanometer
Model	X5
Display Digital LCD	Digital LCD
Method of measurement	Oscillometry
Weight	About 245 g
Dimensions	Appr. 135 (L) X 100 (W) X 68 mm (H)
Cuff	Appr. 150 mm (W) x 450 mm (L)
Cuff Measurement perimeter of the upper arm	220mm~420mm
Material of cuff	Nylon polyester
Measuring pressure range	Blood pressure: 0~295 mmHg, 0~40 kPa
Measuring pressure accuracy	Pressure within ±3 mmHg (±0.4 kPa), or 2% of the reading, whichever is greater
Pulse range	(40~160) times per minute
Pulse accuracy	Within ±5% of the reading
Pressure measurement	Resistance-type pressure sensor
Working mode	Continuous operation
Number of memory	2 people × 99 groups
Air-blow device	Blow-off valve

1.Product Picture

7

8


air tube will be at the side of your elbow. Be careful not to rest your arm on the air tube.

- 1)Sit in a comfortable chair with your back and arms supported, and maintain the right posture; do not force with your fingers, and keep them stretched.
2)Keep your feet flat and your legs uncrossed.
3)Both arms should be resting comfortably on a table.
4)Place your monitor close enough to touch with your bent elbows.
5)The arm cuff should be placed on your arm at the same level as your heart for a blood pressure measurement.

- 6)Please measure while the body is relaxed, and keep it in peace for 4~5 minutes before measurement. Otherwise, the measurement results will be affected. At the same time, make blood pressure measurement with the same arm and the same posture every day.
7)Judge the measurement result based on the guidance of a doctor to help you manage your health.

- 8)When the person being measured suffers from arrhythmia (i.e. PAC, PVC, and atrial fibrillation), the measurement results may be inaccurate.
3.How to Sit Correctly
The correct posture is also extremely important to measurement. To take a measurement, you need to be relaxed and comfortably seated at comfortable room temperature.

- Notes:
• When you take a measurement on the right arm, the

13

6)Avoid any electromagnetic interference or noise interference during measurement.

VIII.Method to measure the blood pressure

- 1)Start measuring after you have worn the cuff by pressing the "Start/Stop" key. After that, start the measurement.

- 2)Since the monitor also checks the pulse during the inflation process, do not move any part of your body during measurement.

Initial display During inflation During measuring

3)When the monitor is slowly vented, the value on the display keeps decreasing and the heartbeat icon (heart shaped) keeps flashing until the measurement results appear.

During inflation During venting During measuring

Reading of measured values and make repeated measurements. After the automatic measurement is finished, the display screen will automatically display the systolic pressure and diastolic pressure of the person under measurement in mmHg, and simultaneously display the measured pulse rate (in times/minute).

14

15

II.Range of application

1.Intended Use

Electronic Sphygmomanometers measures the diastolic pressure and systolic pressure and pulse rate of people over 12 years old in home environments by using a non-invasive oscillometric technique with a single upper arm cuff.

The measurement position is at human being's upper arm.

All values can be read out in one LCD panel. This device is specially designed for family use. It is recommended for people over 12 years old with upper arm circumference between 22cm~42cm.

Scope of application of product
This product is used for measurement of blood pressure of people over 12 years old and not for newborns, not for pregnant, including pre-eclamptic patients;

Intended Users

Adult, except contraindications

Patient populations

People over 12 years, except contraindications

Parts of the body/tissues, duration

This Non-invasive Electronic Sphygmomanometers, the cuff will be worn in upper arm by user during measuring.

Single Use/Reusable
This product is an active and reusable device.

2.Contraindications

1)The product shall not be used for measuring for under 12 years old.

2)Those who suffer from severe arrhythmia shall not use this product;

3)Do not use this monitor on the injured arm or the arm under medical treatment;

4)Do not apply the arm cuff on the arm while being on an intravenous drip or blood transfusion.

III.Important Safety Information

1.Warnings

⚠ To measure the blood pressure accurately and reliably, followings shall be noted when measuring:

- 1)Don't disassemble the monitor or arm cuff. This may cause a serious problem.

- 2)Do not apply the arm cuff on the arm while being on an intravenous drip or blood transfusion.

- 3)Consult your physician before using this monitor on the arm with an arterio-venous (A-V) shunt.

- 4)Don't use the Cuff on the arm on the side of a mastectomy or lymph node clearance.

- 5)Do not use this monitor with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.

- 6)Measurements may be impaired if the device is used close to televisions, microwave ovens, cellular telephones, X-ray or other device with strong electrical fields.

- 7)Don't use this monitor on the injured arm or the arm under medical treatment, as this can cause further injury.

- 8)Do not use in a location with moisture, or a location where water may splash on this monitor. This may damage this monitor.

- 9)Use only the approved arm cuff for this monitor. Use of other arm cuffs may result in incorrect measurement results.

- 10)Do not fold or kink the air tube while taking a measurement. This may cause harmful injury by interrupting blood flow.

- 11)Do not measure while standing, walking, or your body is under pressure.

- 12)Do not smoke after smoking, drinking alcohol, or coffee (black tea).

- 13)Do not measure after exercise or bathing.

- 14)Do not speak or move your body during measurement.

- 2.Notes**

- 1)Avoid eating, smoking and strenuous exercise for 30 minutes prior to blood pressure measurement.
2)Since the blood pressure measurement posture for measurement of blood pressure, and the Cuff and the heart should be on the same level. Therefore, the blood pressure value measured at the front of the counter may be slightly different from that is measured correctly.

- 3)Any blood pressure measurement is affected by the posture and physical condition of the subject. Please be quiet while measuring, and do not move your arms or body.

- 4)Maintain the same state and same posture for each time of measurement.

- 5)The monitor of this model should be tied on the left arm.

4

6)When continuous measurement is made on the same person, loosen the Cuff to have the arm rested for at least 3 minutes before next measurement.

7)The Cuff provided by this product is suitable for arms of people over 12 years old.

8)Avoid electromagnetic interference above 10V/m during measurement. The possible sources of interference may be from TV, mobile phone, microwave oven, X-ray, etc.

9)Do not measure the blood pressure of newborn babies using this product.

10)Do not have this product used for purposes other than blood pressure measurement.

11)Remove the batteries when the product is not in use for a long time to prevent the battery from leakage.

12)When the bladder of cuff is over-pressurized constantly, the pain, numbness and even congestion may occur to the arm. Do not pressurize when it exceeds 300 mmHg, and do not maintain the cuff pressure higher than 15mmHg for more than 3 minutes.

13)Arrhythmia such as PAC, PVC, and atrial fibrillation during the measurement process may lead to the detection of irregular pulse rate. The electronic monitor can still work normally, but the measured blood pressure value may be inaccurate. When irregularities occur frequently, please consult a doctor.

14)If the monitor is used and stored beyond the temperature and humidity range specified in this Manual, its safety or performance may be affected or become invalid, or even be totally damaged.

15)When replacing the Cuff, use the Cuff supplied by the manufacturer. If the original part is replaced with a Cuff not supplied by the manufacturer, measurement errors may occur.

16)Contact the manufacturer for static pressure test.

17)Do not share the Cuff with others to avoid cross-contamination.

5


VII.Preparation

Open the Electronic Sphygmomanometers. The device case should be uniformly colored with no obvious defects such as scars and scratches.

Please keep it properly when carrying and storing.

1.To load the batteries

- 1)Open the Battery Case cover on the side of the unit and put the 4 Batteries (AAA) correctly.
2)Put the Battery Case cover back properly.

Notes:

When the depleted battery symbol () appears on the display, turn the monitor off and remove all the batteries. Replace with 4 new batteries at the same time. Low voltage alkaline batteries should be discarded.

The measurement values continue to be stored even after the batteries are replaced.

The supplied batteries may have a shorter life. Disposal of used batteries should be carried out in accordance with the national/local regulations for the disposal of batteries.

2.To tie the Cuff

- 1)The standard accessory is the cuff (M), which is suitable for the subject with arm circumference of 22cm ~ 42cm.
2)Insert the hose top plug into the opening located on the side of the Main Unit. (Only the original standard Cuff can be used for this Main Unit.)

12

X.I Requirements of Storage and maintenance

1.Avoid falling or strong collisions.

2.Avoid high temperatures or exposure to direct sunlight.

3)Don't twist or bend the Cuff.

4)Do not press "START/STOP" button if the cuff is not wrapped around your upper

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